

# Development of type 2 and type 3 iodothyronine deiodinase knockout *Xenopus tropicalis* using CRISPR/Cas12a gene editing

Sally Mayasich<sup>1,2\*</sup>, Philip DeGoey<sup>2</sup>, Jonathan Haselman<sup>2</sup>, Sigmund Degitz<sup>2</sup>


(1) ORISE Post-doctoral fellow, (2) U.S. Environmental Protection Agency, Office of Research and Development, Center for Computational Toxicology and Exposure, Great Lakes Toxicology and Ecology Division, Duluth, MN, USA

*Disclaimer: The views expressed in this presentation are those of the authors and do not necessarily reflect the views or policies of the U.S. Environmental Protection Agency. Mention of trade names or commercial products does not constitute endorsement or recommendation for use.*

*\*Current affiliation: University of Wisconsin-Madison Cooperative Post-doctoral fellow at U.S. EPA GLTED*


# Background

- Iodothyronine deiodinase enzymes (dio1, 2, and 3) catalyze tissue-localized metabolic activation/deactivation of thyroid hormones (THs)
- We targeted the dio2 (TH activating) and dio3 (TH inactivating) enzymes in this study



# Postulated pathway for amphibian deiodinase gene disruption

**Purpose of genetic modifications:** this work supports the EPA's chemical safety research mission by using an alternative approach for linking molecular mechanisms of thyroid toxicity to apical outcomes relevant to ecological risk assessment



## Overview


- We utilized a CRISPR/Cas12a system to induce DNA cleavage in the *dio* genes at locations near where the catalytic sites would occur in the translated proteins
  - The Cas12a Protospacer Adjacent Motif (PAM) sites gave us better cutting locations for these genes than those offered by the Cas9 system
  - We also got slightly better efficiency with Cas12a than Cas9 in tests to knock out the tyrosinase gene
  - Cas12a is more convenient; requires only single guide RNA (no tracrRNA)
- Single-cell-stage zygotes were injected with either ***dio2*** or ***dio3*** gene-specific guide-RNA/Cas12a ribonucleo-protein (RNP) complexes

## CRISPR/Cas12a (LbCpf1) crRNA design and microinjection methods

- Guide RNA design was aided using CRISPRscan and CRISPR RGEN tools Cas-designer
- The CRISPR guide RNA oligos (crRNA) were synthesized by Integrated DNA Technologies and tested through in vitro cleavage with EnGen Lba Cas12a (LbCpf1) DNA endonuclease (New England BioLabs)
- *X. tropicalis* adults were induced to spawn by hCG injection; eggs were treated with L-cysteine rinse; single-cell zygotes were sorted into 1/9x Modified Ringers with 3% Ficoll
- Zygotes were injected with ~2 nl (20 fmol) of either the dio2 or dio3 RNP complex in Cpf1 working buffer (20mM Hepes pH 7.5, +KCl+MgCl<sub>2</sub>+TCEP+glycerol) Moreno-Mateos et al. (2017)

<https://doi.org/10.1038/s41467-017-01836-2>

## crRNA-LbCpf1 Ribonucleoprotein complex for dio3





**TGA/UGA** = Codon for critical selenocysteine amino acid in deiodinase enzyme active site.

**PAM** = Protospacer adjacent motif



## Microinjection of RNP complexes into *X. tropicalis* zygotes




## Culture of injected *X. tropicalis* tadpoles

Following 24-30 hours of incubation at 25°C in 1/9x modified Ringers, embryos were transferred to 4.5 L flow-through tanks delivering Lake Superior Water @ ~25 mL/min.





# DNA extraction and polymerase chain reaction (PCR) genotype screening



Agilent Bioanalyzer capillary electrophoresis. Crisprants (animals with disrupted genes) show multiple bands due to heteroduplex formation during PCR.

## Results

- PCR genotyping at 3 weeks post-injection revealed putative mutations at a rate of >60% in dio2-RNP-injected tadpoles, and 50% in dio3-RNP-injected tadpoles
- Apical outcomes included asynchronous metamorphosis in dio2 and high mortality in dio3 crispants


## Apical outcomes of deiodinase gene disruption

Injected, unaffected **WT** at 45 days  
Nieuwkoop and Faber (NF) stage 66




***dio2*** crispant at 44 days  
Asynchronous metamorphosis  
(differential limbs)

***dio3*** crispant at 26 days  
Precocious metamorphosis/early  
tail and gill resorption




## Metamorphosis timing in *dio2* crispant and WT phenotypes


### *X. tropicalis* *dio2* CRISPR Days to Stage 58



### *X. tropicalis* *dio2* CRISPR Days to Stage 62




## Precocious metamorphosis leads to high mortality in *dio3* crispants






## Cross breeding conducted



## Example *dio2* crispant cross F<sub>1</sub> genotyping results and interpretation hypothesis



### Summary:

- Genotyping results show indel mutations occurred in CRISPR/Cas12a *dio* gene-edited frogs;
- Crispant F0 apical outcomes indicate putative loss of *dio* activity leading to altered metamorphosis.

### Next steps:

- Genotyping of *dio3* crosses and sequencing;
- Evaluation of F1 apical outcomes;
- Multiple *dio* knockouts?

## Acknowledgements

### **EPA/GLTED Thyroid team**

Jessica Christensen

Jeff Denny

Thomas Lerdall

Joe O'Flanagan

Jennifer Olker

### **Badger Technical Services**

Stephanie Wisniewski

This work was funded wholly by the U.S. Environmental Protection Agency.

*Disclaimer: The views expressed in this presentation are those of the authors and do not necessarily reflect the views or policies of the U.S. Environmental Protection Agency. Mention of trade names or commercial products does not constitute endorsement or recommendation for use.*