

High-throughput phenotypic profiling within the NAMs-based, tiered hazard evaluation strategy at the United States Environmental Protection Agency

Johanna Nyffeler


ORISE postdoctoral grantee at USEPA Center for Computational Toxicology and Exposure (CCTE)

ORCID 0000-0002-6155-9743

Nyffeler.Johanna@epa.gov

The views expressed in this presentation are those of the author and do not necessarily reflect the views or policies of the U.S. EPA.

Tiered Hazard Evaluation Strategy based on New Approach Methods (NAMs)




High-Throughput Phenotypic Profiling (HTPP)

- labeling of various cell organelles with fluorescent probes in *in vitro* cultures
- assessing a large variety of morphological features on individual cells via imaging

'Cell Painting' assay

Gustafsdottir *et al.* 2013,
PMID: 24312513

Bray *et al.* 2016,
PMID: 27560178



Fluorescent labels

DNA: H-33342

RNA: SYTO14

ER: Concanavalin A-488

Actin: Phalloidin-568



Golgi + Membrane: wheat
germ agglutinin (WGA) -555

Mitochondria: MitoTracker

Nyffeler *et al.* 2020




HTPP: Two Applications




← profile

for each chemical x concentration



Application 1

concentration-response modelling



Potency estimation:
→ *in vitro* point-of-departure (POD)

Application 2

Chemical A

0	1.80	0	9.69	5.73	0	6.47	-12.84	0	0
---	------	---	------	------	---	------	--------	---	---

Biological similarity

Chemical B

0	0	0	10.00	6.00	1.60	6.47	-15.00	0	0
---	---	---	-------	------	------	------	--------	---	---

Compare profiles with annotated reference chemicals
→ putative mechanisms

- Nyffeler *et al.* (2020) Toxicol Appl Pharmacol. PMID: 31899216
- Willis *et al.* (2020). SLAS Discov. PMID: 32546035
- Nyffeler *et al.* (2021). SLAS Discov. PMID: 32862757


work in progress




HTPP Screening Results

Application 1:

Potency estimation




Active chemicals:



- ⇒ ~ 40% of chemicals were active
- ⇒ Most activity is > 10 μM

Comparison with ToxCast screening results:




⇒ **Less potent than ToxCast POD**



Comparison to *in vivo* Effect Values

Application 1:

Potency estimation



- 303 chemicals were active and had pharmacokinetic (PK) information





⇒ **78% of HTPP AED are within 2 orders of magnitude of the *in vivo* POD**

POD: point-of-departure
AED: administered equivalent dose


Comparison to Exposure Estimates

Application 1:

Potency estimation



Potential for humans
to be exposed to
bioactive concentrations




- ⇒ for 49% of chemicals, predicted exposure is > 1000x lower than estimated bioactivity
- ⇒ for a small set of chemicals, the BER was negative, indicating a potential for humans to be exposed to bioactive concentrations of these chemicals



Preliminary results. Do not cite or quote.

Feature Selection & Profile Comparison

Application 2: Mechanistic prediction




Feature Selection

1300 features




317 features

Profile Comparison



Similar Mechanism → Similar Phenotype

Application 2: Mechanistic prediction




- 50 chemicals were annotated as targeting a nuclear receptor


target

AR
ESR
GR
PGR
PPAR
PXR
RAR
VDR

Biological similarity in HTPP
50 chemicals



Gene expression in U-2 OS




⇒ Agonists of the GR and of RAR/RXR display characteristic profiles

Preliminary results. Do not cite or quote.


Application: Find Retinoid-like Chemicals

Application 2:

Mechanistic prediction




1. Compare profiles to 5 known retinoids:



→ 10 candidate chemicals


2. Repeat HTPP experiments:



→ 4/10 chemicals confirmed

3. Compare to ToxCast results:

known retinoids




→ These 4 chemicals had activity in ToxCast assays targeting RAR/RXR


⇒ HTPP has the potential to identify environmental chemicals with specific activities

Specific vs Non-specific Phenotypes

Application 2: Mechanistic prediction




small groups of
chemicals with a
very **specific**
phenotype



chemicals (n = 470)

chemicals (n = 470)

large groups of
chemicals with
similar phenotype
→ **non-specific**
toxicity
mechanisms?




Application: Grouping of Conazoles


Application 2: Mechanistic prediction

- group of fungicides
- disturb ergosterol synthesis via CYP51 and CYP61 (target absent in mammals)

biological similarity




structural similarity (based on ToxPrints)



- ⇒ most conazoles are phenotypically similar
- ⇒ Diniconazole is phenotypically different from the other active conazoles

Conclusions




Application 1: Potency estimation

- HTPP can be used to derive *in vitro* potency estimates
- *in vitro* potency estimates often comparable or more conservative than *in vivo* PODs
- used for Bioactivity-Exposure-Ratio (BER) analysis

Application 2: Mechanistic prediction

- Identification of chemicals with specific mechanisms e.g., chemicals with retinoid-like activity
- Biological grouping of structurally related chemicals e.g., conazoles



Acknowledgements

Office of Research and Development (ORD) Center for Computational Toxicology and Exposure (CCTE)



Harrill Lab team

- **Joshua Harrill**
- **Clinton Willis**
- **Felix Harris**
- Rick Brockway
- Megan Culbreth
- Dan Hallinger
- Terri Fairley

Data analysis

- Daniel Chang
- Kathy Coutros
- Logan Everett
- Katie Paul Friedman
- Derik Haggard
- Richard Judson
- Ryan Lougee
- Grace Patlewicz
- Ann Richard
- Woody Setzer
- Imran Shah
- John Wambaugh

