



Methylene Blue and Developing Zebrafish – Effects on Mortality, Morphology and Behavior

Joan M. Hedge

Office of Research and Development


Center for Computational Toxicology and Exposure/ Biomolecular and Computational Toxicology Division/Advanced Experimental Toxicology Models Branch

Research Triangle Park, North Carolina



This presentation does not necessarily reflect official United States Environmental Protection Agency policy and mention of trade names or commercial products does not constitute endorsement or recommendation for use.



Does your embryo media with methylene blue added look like a darker blue than this?





If so, your zebrafish may be in trouble!



Methylene Blue




- “First fully synthetic drug used in medicine”
 - Humans - treatment of methemoglobinemia, marking for polyps/tumor removal and to show lymphatic drainage
 - Humans – previously used for cyanide poisoning, urinary tract infections, malaria, bipolar disorder, Alzheimer’s disease, RNA viruses, and marking the amniotic sac
 - Mice - shown to cause both axial skeleton and neural tube defects
 - Aquaculture/research and tropical fish hobbyists - treatment for fungal infections and to protect newly fertilized fish embryos from infection
- 




- Use on zebrafish embryos appears to be widespread
- Reported concentrations vary - 0.03 µM to over 3100 µM (husbandry manuals, experimental protocols & research papers)
- More commonly used concentrations on zebrafish embryos are 10 µM or less
- No systematic study of toxicological effects on developing zebrafish / no standardization for use

Experimental Protocol


Pre-assay




During assay



Post-assay




Embryo Mortality




Control n=36

Doses 10-250 µM, n=24

- Embryo survival was not impacted by methylene blue in concentrations up to 250 µM




Dead, 6 dpf



Alive, 6 dpf

Fisher's Exact Test on Incidence Data, $P > 0.05$

Embryo Hatching



Control n=36

Doses 10-250 µM, n=24

Fisher's Exact Test on Incidence Data, $P < 0.0055$ (125µM); $P < 0.00001$ (250 µM)

- Reduction in percent of embryos hatched with concentrations of $\geq 125 \mu\text{M}$ methylene blue




Unhatched, 6 dpf



Hatched, 6 dpf

Percent Normal Embryos



Control n=36

Doses 10-250 µM, n=24

Fisher's Exact Test on Incidence Data, $P < 0.00001$

- Removing dead, unhatched, abnormal (swim bladder uninflation, etc.) embryos, % “normal” was calculated
- Reduction in the percent of “normal” embryos with 250 µM methylene blue
- No “normal” embryos found at 250 µM methylene blue




Uninflated swim bladder, 6 dpf



Normal, 6 dpf


Body Length



Control n=36


Doses 10-250 µM, n=10-24

- No dose response relationship between methylene blue and body length.



ANOVA, P-Value = 0.0571


Maximum Tail Curvature



Control n=36


Doses 10-250 μM, n=10-24

- Maximum tail curvature increased with 250 μM methylene blue



ANOVA, P-Value < 0.0001, followed by Fisher's PLSD


Swim Bladder Size



Control n=36


Doses 10-250 µM, n=10-24

- Swim bladder size reduced in larvae with 125 µM or > methylene blue




ANOVA, P-Value < 0.0001, followed by Fisher's PLSD

Eye Size



Control n=36 Doses 10-250 µM, n=10-24 ANOVA, P-Value < 0.0001, followed by Fisher's PLSD

- Eye size = most sensitive endpoint assessed
- 63 µM or > methylene blue caused a significant reduction in eye size



Effects of Methylene Blue

	Locomotor Behavior	Mortality	Hatching	% Normal	Body Length	Tail Curvature	Swim Bladder Size	Eye Size
less than 10 μM	-	-	-	-	ND	ND	ND	ND
10 μM	ND	-	-	-	-	-	-	-
31 μM	ND	-	-	-	-	-	-	-
63 μM	ND	-	-	-	-	-	-	+
94 μM	ND	-	-	-	-	-	-	+
125 μM	ND	-	+	-	-	-	+	+
250 μM	ND	-	+	+	-	+	+	+

0 μM

10 μM

31 μM

63 μM


94 μM

125 μM

250 μM

Conclusions

- Including our previous work, methylene blue appears safe for use with zebrafish embryos up to **10 μM** which falls within the most commonly used concentrations.
- Concentrations above **31 μM** methylene blue are not recommended.
- Use of methylene blue when the exact concentration is not known is not safe for use with zebrafish embryos.



THANK YOU.....



Contributors to this research

- Stephanie Padilla
- Erik Sanders
- Jeanene Olin
- Bridgett Hill
- Deborah Hunter
- Morgan Lowery
- Kimberly Jarema
- Alan Tennant
- Susan Farmer

Special thanks to the animal care staff

- Donald Holman
- Guillermo Orozco
- Femi Yerumo
- James Cicala
- Mark King
- Leslie Jarrell
- Kim Howell
- Ned Collins