




www.epa.gov

Background

- The Human Pluripotent Stem Cell Test (hPST) utilizes stem cell endoderm differentiation to identify teratogens that perturb SOX17 expression (1).
- To adopt the principles of the hPST, we utilized the RUES2-GLR (germ layer reporter) cell line engineered with fluorescent reporters for biomarkers representing all three gastrulation lineages: endoderm, ectoderm and mesoderm (2).
- This RUES2-GLR assay enables rapid 384-well high-throughput screening and quantitative analysis of multi-lineage differentiation, demonstrating potential to identify and prioritize teratogens based on affected gastrulation lineages.

Method

- Human pluripotent stem cell line (RUES2-GLR) containing a fluorescent biomarker, SOX17–tdTomato (endoderm), was used to assess perturbations during endoderm differentiation.



- High-content image analysis was used to quantitatively determine the percentage of SOX17+ cells and total cell count (viability).

References

- 1) Sei Kameoka, Joshua Babiarz, Kyle Kolaja, Eric Chiao, A High-Throughput Screen for Teratogens Using Human Pluripotent Stem Cells, Toxicological Sciences, Volume 137, Issue 1, January 2014, Pages 76–90.
- 2) Martyn, I., Kanno, T.Y., Ruza, A. et al. Self-organization of a human organizer by combined Wnt and Nodal signalling. Nature 558, 132–135 (2018).

This poster does not necessarily reflect EPA policy.

U.S. Environmental Protection Agency
Office of Research and Development


Adaptation of a Human Pluripotent Stem Cell Assay for Developmental Toxicity Screening

John Gamble^{1,2}, Kristen Hopperstad¹, Chad Deisenroth¹

¹Center for Computational Toxicology and Exposure, US EPA Research Triangle Park, NC, USA


²Oak Ridge Institute for Science and Education, Oak Ridge, TN, USA

Email: gamble.john@epa.gov



Inhibition of SOX17 expression during endoderm induction is more predictive of a developmental toxicant than cytotoxicity.

Evaluation of assay predictivity using a 66 chemical training set with known developmental effects.



Assay Performance

- rZ' : 0.78
- Intra rCV: 6.4%
- Inter rCV: 5.3%
- Sensitivity: 79%
- Specificity: 65%
- Bal. accuracy: 72%